

SCIENTIFIC COOLING		17th Feb 2025	
Pharmacy fridges compliant with DIN 58345		RRP excl. vat	RRP incl. vat
HMTvh 1501 (Replacement for MKUv 1610)	Fridge Mediline +5°C, automatic defrost, dynamic cooling, electronic controller, external digital temp display, power failure alarm immediately upon power failure for 12hours, 3 adjustable shelves, self-closing door, lockable, white finish, solid door, size 597mm(w) x 605mm (d) 819mm capacity 102l nett	R 55 643,48	R 63 990,00
HMTvh 1511 (Replacement for MKUv 1613)	Fridge Mediline +5°C, automatic defrost, dynamic cooling, electronic controller, external digital temp display, power failure alarm immediately upon power failure for 12hours, 3 adjustable shelves, self-closing door, lockable, white finish, glass door, size 597mm(w) x 605mm (d) 819mm (h), capacity 102l nett	R 58 252,17	R 66 990,00
HMfvh 4001 (Replacement for MKv 3910)	Fridge Mediline +5°C, automatic defrost, dynamic cooling, electronic controller, external digital temp display, power failure alarm immediately upon power failure for 12hours, 5 adjustable shelves, self-closing door, lockable, white finish, solid door, size 597mm(w) x 654mm (d) 1660mm (h), capacity 295l nett	R 61 730,43	R 70 990,00
HMfvh 4011 (Replacement for MKv 3913)	Fridge Mediline +5°C, automatic defrost, dynamic cooling, electronic controller, external digital temp display, power failure alarm immediately upon power failure for 12hours, LED interior light with separate switch, 5 adjustable shelves, self-closing door, lockable, white finish, glass door, size 597mm(w) x 654mm (d) 1884mm (h), capacity 297l nett	R 65 208,70	R 74 990,00
Laboratory fridges and freezers with Profi electronic controller		RRP excl. vat	RRP incl. vat
SRPvg 1412 (Replacement for LKPv 1423)	Fridge Mediline 0°C to +16°C, automatic defrost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, LED interior light with separate switch, 8 adjustable shelves, self-closing door, lockable, white finish, glass door, castors, size 1427mm(w) x 867mm (d) 2135mm (h), capacity 1447l gross	R 156 513,04	R 179 990,00
SRPvg 1402 (Replacement for LKPv 1420)	Fridge Mediline -2°C to +16°C, automatic defrost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, 8 adjustable shelves, self-closing door, lockable, white finish, solid door, castors, size 1427mm(w) x 867mm (d) 2135mm (h), capacity 1391l gross	R 147 817,39	R 169 990,00
SRPvg 8401 (Replacement for LKPv 8420)	Fridge Mediline -2°C to +16°C, automatic defrost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, 4 adjustable shelves, self-closing door, lockable, white finish, solid door, castors, size 787mm(w) x 1017mm (d) 2135mm (h), capacity 871l gross	R 104 339,13	R 119 990,00
SRPvg 6511 (Replacement for LKPv 6523)	Fridge Mediline 0°C to +16°C, automatic defrost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, LED lighting, 4 adjustable shelves, self-closing door, lockable, white finish, glass door, castors, size 697mm(w) x 867mm (d) 2135mm (h), capacity 642l gross	R 104 339,13	R 119 990,00
SRPvg 6501 (Replacement for LKPv 6520)	Fridge Mediline -2°C to +16°C, automatic defrost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, 4 adjustable shelves, self-closing door, lockable, white finish, solid door, castors, size 697mm(w) x 867mm (d) 2135mm (h), capacity 614l gross	R 101 730,43	R 116 990,00
SFPvg 1402 (Replacement for LGPv 1420)	Freezer Mediline -9°C to -30°C, automatic defrost, no frost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, LED interior light with separate switch, 8 adjustable shelves, self-closing door, lockable, white finish, solid door, castors, size 1427mm(w) x 867mm (d) 2135mm (h), capacity 1391l gross	R 147 817,39	R 169 990,00
SFPvg 8401 (Replacement LGPv 8420)	Freezer Mediline -9°C to -30°C, automatic defrost, no frost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, LED ceiling lighting, 4 adjustable shelves, self-closing door, lockable, white finish, solid door, castors, size 787mm(w) x 1017mm (d) 2135mm (h), capacity 871l gross	R 99 991,30	R 114 990,00
SFPvg 6501 (Replacement for LGPv 6520)	Freezer Mediline -9°C to -30°C, automatic defrost, no frost, dynamic cooling, Profi electronic controller with data logger and text display, external digital temp display, power failure alarm immediately upon power failure for 72hours, LED ceiling lighting 4 adjustable shelves, self-closing door, lockable, white finish, solid door, castors, size 697mm(w) x 867mm (d) 2135mm (h), capacity 614l gross	R 99 991,30	R 114 990,00
Laboratory fridges and freezers with Comfort electronic controller		RRP excl. vat	RRP incl. vat
SRFvg 4011 (Replacement model for LKV 3913)	Fridge Mediline +3°C to +16°C, automatic defrost, dynamic cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 5 adjustable shelves, self-closing door, lockable, white finish, LED lighting, glass door, size 597mm(w) x 654mm (d) 1884mm (h), capacity gross/net 420/297l	R 54 773,91	R 62 990,00
SRFvg 4001 (Replacement model for LKV 3910)	Fridge Mediline +3°C to +16°C, automatic defrost, dynamic cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 5 adjustable shelves, self-closing door, lockable, white finish, solid door, size 597mm(w) x 654mm (d) 1884mm (h), capacity gross/net 394/298l	R 52 165,22	R 59 990,00
SRTfg 1501 (Replacement model for LKUxv 1610)	Fridge Mediline +3°C to +16°C, Spark-free interior, automatic defrost, dynamic cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 3 adjustable shelves, self-closing door, lockable, white finish, solid door, size 601mm(w) x 618mm (d) 819mm (h), capacity 102l nett	R 47 817,39	R 54 990,00
Added Storage fees			
	Month 1-3 No Charge		
	Month 4-6 Per Unit / Per Month	R 65,22	R 75,00
	Month 7-12 Per Unit / Per Month	R 86,96	R 100,00
	Month 13 Onwards	R 130,43	R 150,00

"NB Due to the instability of our National Electricity Supplier, Liebherr strongly recommends that all appliances are fitted with a surge protector, as instances of this nature, are not considered a factory fault."

SCIENTIFIC COOLING		17th Feb 2025	
SRTvg 1511 (Replacement model for LKUV 1613)	Fridge Mediline +3°C to +16°C, automatic defrost, dynamic cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, LED ceiling light, 3 adjustable shelves, self-closing door, lockable, white finish, glass door, size 597mm(w) x 605mm (d) 819mm (h), capacity 102l nett	R 50 426,09	R 57 990,00
Laboratory fridges and freezers with Comfort electronic controller		RRP excl. vat	RRP incl. vat
SRTvg 1501 (Replacement model for LKUV 1610)	Fridge Mediline+3°C to +16°C, automatic defrost, dynamic cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 3 adjustable shelves, self-closing door, lockable, white finish, solid door, size 597mm(w) x 605mm (d) 819mm (h), capacity 102l nett	R 46 078,26	R 52 990,00
Laboratory fridge / freezers with Comfort electronic controller		RRP excl. vat	RRP incl. vat
SCFvh 4002 (Replacement for LCV 4010)	Fridge Freezer Mediline +3°C to +16°C / -9°C to -30°C, automatic defrost, no frost, forced air/static cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 4 adjustable shelves, self-closing door, lockable, LED Interior, white finish, solid doors, size 597mm(w) x 654mm (d) 2044mm (h), capacity 272l / 206l gross/nett	R 74 773,91	R 85 990,00
Spark-free Laboratory fridges and freezers with Comfort electronic controller		RRP excl. vat	RRP incl. vat
SCFFg 4002 (Replacement for Lcexv 4010)	Fridge Freezer Mediline+3°C to +16°C / -9°C to -30°C, Spark-free interior, smart frost, automatic defrost, forced-air / static cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 4 adjustable shelves, 3 drawers, self-closing door, lockable, steel / white finish, solid door, size 597mm(w) x 654mm (d) 2044mm (h), capacity fridge 199l / freezer 62l net	R 63 469,57	R 72 990,00
SRFfg 4001 (Replacement model for Lkexv 3910)	Fridge Mediline +3°C to +16°C, Spark-free interior, automatic defrost, dynamic cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 5 adjustable shelves, self-closing door, lockable, white finish, solid door, size 597mm(w) x 654mm (d) 1884mm (h), capacity gross/net 394/298l	R 49 556,52	R 56 990,00
SFFfg 4001 (Replacement model for LGEX 3410)	Freezer Mediline -9°C to -30°C, Spark-free interior, manual defrost, smart frost, static cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 5 adjustable shelves, lockable, white finish, solid door, size 597mm(w) x 654mm (d) 1884mm (h) capacity gross/net 316l / 242l	R 57 382,61	R 65 990,00
SFTfg 1501 (Replacement for LGUex 5501)	Freezer Mediline -9°C to -30°C, Spark-free interior, manual defrost, smart frost, static cooling, Comfort electronic controller, external digital temp display, power failure alarm when power returns, 2 shelves, lockable, white finish, solid door, size 601mm(w) x 618mm (d) 820mm (h), capacity 73l nett	R 44 339,13	R 50 990,00
Spark-free Laboratory fridges and freezers with mechanical control		RRP excl. vat	RRP incl. vat
SRFfg 5501 (Replacement model for Lkexv 5400)	Fridge Mediline +3°C to 16°C, Spark-free interior, automatic defrost, dynamic cooling, analogue controller, external digital temp display, power failure alarm when power returns, 5 adjustable shelves, lockable, white finish, solid door, size 747mm(w) x 769mm (d) 1684mm (h) capacity gross/net 558/441l	R 52 165,22	R 59 990,00
SRFfg 3501 (Replacement model for Lkexv 3600)	Fridge Mediline +3°C to 16°C, Spark-free interior, automatic defrost, dynamic cooling, analogue controller, external digital temp display, power failure alarm when power returns, 4 adjustable shelves, lockable, white finish, solid door, size 597mm(w) x 654mm (d) 1684mm (h), capacity gross/net 344/261l	R 47 817,39	R 54 990,00
"NB Due to the instability of our National Electricity Supplier, Liebherr strongly recommends that all appliances are fitted with a surge protector, as instances of this nature, are not considered a factory fault."			
Added Storage fees	Month 1-3 No Charge		
	Month 4-6 Per Unit / Per Month	R 65,22	R 75,00
	Month 7-12 Per Unit / Per Month	R 86,96	R 100,00
	Month 13 Onwards	R 130,43	R 150,00