

If you cannot find your appliance in the list below, please visit [our warrantyplus page](#) and enter your serial number, or email us at LGB.websales@liebherr.com

American style fridge-freezers –

15% discount when purchasing WarrantyPlus for both appliances at the same time

Freestanding

SBS Model	Sides	+3 years	+8 years
SBSes 8496	SKBes 4380	£ 150	£ 300
	SWTNes 4285	£ 220	£ 440
XCCsd 5250	SCNsdd 5253/SCNsdd 525i	£ 185	£ 370
	SCNsdd 5253/SCNsdd 525i	£ 185	£ 370
XRCsd 5255	SBNsdd 5264/SBNsdd 526i	£ 185	£ 370
	SRsdd 5250/SRsdc 525i	£ 125	£ 250
XRF 5220	SFNe 5227/SFNd 5227	£ 105	£ 210
	SRe 5220/SRd 5220	£ 105	£ 210
XRFbd 5220	SFNbde 5227/SFNbdd 5227	£ 105	£ 210
	SRbde 5220/SRbdd 5220	£ 105	£ 210
XRFbs 5295	SFNbsd 529i	£ 150	£ 300
	SRBbsd 529i/SRBbsc 529i	£ 150	£ 300
XRFsd 5220	SFNsde 5227	£ 105	£ 210
	SRsde 5220	£ 105	£ 210
XRFsf 5220	SFNsfe 5227/SFNsfd 5227	£ 105	£ 210
	SRsfe 5220/SRsfd 520	£ 105	£ 210
XRFsf 5225	SFNsfe 5227/SFNsfd 5227	£ 105	£ 210
	SRBsfe 5220/SRBsfc 5220	£ 105	£ 210
XRFsf 5240	SFNsfe 5247/SFNsfd 5247	£ 105	£ 210
	SRsfe 5220/SRsfd 520	£ 105	£ 210
XRFsf 5245	SFNsfe 5247/SFNsfd 5247	£ 105	£ 210
	SRBsfe 5220/SRBsfc 5220	£ 105	£ 210
XRFsd 5255	SRBsdd 5250	£ 125	£ 250
	SFNsdd 5257	£ 105	£ 210
XRFst 5295	SFNstd 529i	£ 150	£ 300
	SRBstd 529i/SRBstc 59i	£ 150	£ 300

Fridge-freezers

Freestanding

Model	+ 3 years	+ 8 years
CBNa 572i	£160	£320
CBNamy 575i	£185	£370
CBNbda 5223	£160	£320
CBNbda 5723	£160	£320
CBNbda 572i	£160	£320
CBNbdc 573i	£160	£320
CBNbmy 579i	£220	£440
CBNbsa 5753	£185	£370
CBNbsa 575i	£185	£370
CBNbsc 778i	£220	£440
CBNbsd 578i	£220	£440
CBNc 5223	£160	£320
CBNc 5723	£160	£320
CBNcmy 572i	£160	£320
CBNef 5735	£220	£440
CBNes 5775	£220	£440
CBNes 6256	£285	£570
CBNsda 5223	£160	£320
CBNsda 5723	£160	£320
CBNsda 572i	£160	£320
CBNsda 5753	£185	£370
CBNsda 575i	£185	£370
CBNsdb 5753	£185	£370
CBNsdb 575i	£185	£370
CBNsdb 775i	£185	£370
CBNsdc 522i	£160	£320
CBNsdc 573i	£160	£320
CBNsdc 765i	£185	£370
CBNsfc 5223	£160	£320
CBNsfc 522i	£160	£320
CBNsfc 57vi	£160	£320
CBNsfd 5723	£160	£320
CBNstb 579i	£220	£440
CBNstc 579i	£220	£440
CBNste 8872	£285	£570
CNamy 5723	£160	£320
CNbda 5723	£160	£320
CNbdb 5223	£160	£320
CNbdc 573i	£160	£320
CNbdd 5733	£160	£320
CNc 5023	£160	£320

Built-In

Model	+ 3 years	+ 8 years
ECBN 5066	£285	£570
ECBN 6156	£285	£570
ECBN 6256	£285	£570
ICBb 5152	£185	£370
ICBbi 5122	£160	£320
ICBbi 5152	£185	£370
ICBc 5182	£220	£440
ICBci 5182	£220	£440
ICBNci 5153	£185	£370
ICBNci 5183	£220	£440
ICBNcx 5173	£220	£440
ICBNd 5153	£185	£370
ICBNd 5163	£185	£370
ICBNdi 5123	£160	£320
ICBNdi 5163	£185	£370
ICBNdi 5183	£220	£440
ICBNe 5123	£160	£320
ICBNSd 5123	£160	£320
ICBSd 5122	£160	£320
ICc 5123	£160	£320
ICd 5103	£85	£170
ICd 5123	£160	£320
ICe 5103	£85	£170
ICNc 5123	£160	£320
ICNci 5173	£220	£440
ICNd 5123	£160	£320
ICNdi 5153	£185	£370
ICNdi 5173	£220	£440
ICNe 5103	£85	£170
ICNe 5133	£160	£320
ICNf 5103	£85	£170
ICNSf 5103	£85	£170
ICSe 5103	£85	£170
SICNd 5153	£185	£370
SICNdi 5153	£185	£370

Fridge-freezers

Freestanding

Model	+ 3 years	+ 8 years
CNc 5703	£85	£170
CNc 5724	£160	£320
CNcmy 5203	£85	£170
CNcmy 5724	£160	£320
CNd 5203	£85	£170
CNd 5204	£85	£170
CNd 5703	£85	£170
CNd 5704	£85	£170
CNd 5724	£160	£320
CNd 7723	£160	£320
CNsda 5723	£160	£320
CNsdb 5223	£160	£320
CNsdb 5723	£160	£320
CNsdc 5203	£85	£170
CNsdd 775i	£185	£370
CNsdd 776i	£185	£370
CNsfc 5023	£160	£320
CNsfc 573i	£160	£320
CNsfc 574i	£160	£320
CNsfd 5023	£160	£320
CNsfd 5203	£85	£170
CNsfd 5204	£85	£170
CNsfd 5703	£85	£170
CNsfd 5704	£85	£170
CNsfd 7723	£160	£320

Freestanding (continued)

Model	+ 3 years	+ 8 years
CU 2331	£85	£170
CU 2831	£85	£170
CU 3331	£85	£170
CUb 2331	£85	£170
CUb 2331	£85	£170
CUbe 2331	£85	£170
CUe 2331	£85	£170
CUe 2831	£85	£170
CUe 3331	£85	£170
CUef 3331	£85	£170
CUefe 3331	£85	£170
CUel 2331	£85	£170
CUel 2831	£85	£170
CUel 3331	£85	£170
CUele 2331	£85	£170
CUele 2831	£85	£170
CUele 3331	£85	£170
CUfb 2831	£85	£170
CUfbe 2831	£85	£170
CUfr 2331	£85	£170
CUfre 2331	£85	£170
CUno 2831	£85	£170
CUnoe 2831	£85	£170

Top mounted freezer

Model	+ 3 years	+ 8 years
CT 2131	£85	£170
CT 2531	£85	£170
CTe 2131	£85	£170
CTe 2531	£85	£170
CTe 2931	£85	£170
CTel 2131	£85	£170
CTel 2531	£85	£170
CTel 2931	£85	£170
CTele 2131	£85	£170
CTele 2531	£85	£170
CTele 2931	£85	£170

Tall fridges

Freestanding

Model	+ 3 years	+ 8 years
K 2340	£95	£190
K 2834	£95	£190
Ke 2834	£95	£190
Kel 2834	£95	£190
Kele 2834	£95	£190
RBa 4250	£125	£250
RBa 425i	£125	£250
RBbmy 525i	£125	£250
RBbsb 525i	£125	£250
RBbsc 5250	£125	£250
RBc 5220	£105	£210
RBc 525i	£125	£250
RBd 5250	£125	£250
RBsdc 525i	£125	£250
RBsdd 5250	£125	£250
RBsfc 5220	£105	£210
RBsfd 5221	£105	£210
RBstd 528i	£150	£300
Rd 4200	£95	£190
Rd 4600	£95	£190
Rd 5000	£95	£190
Rd 5220	£105	£210
RDsfd 5220	£105	£210
RDsfe 5220	£105	£210
Re 5220	£105	£210
Rsfed 4600	£95	£190
Rsfed 5000	£95	£190
Rsfed 5220	£105	£210
Rsfe 5020	£105	£210
Rsfe 5220	£105	£210

Built-In

Model	+ 3 years	+ 8 years
IRBAc/d 5171	£150	£300
IRBAc/d 5190	£150	£300
IRBb 4170	£150	£300
IRBbi 4050	£125	£250
IRBbi/bx 4170	£150	£300
IRBbi/bx 4570	£150	£300

Built-In

Model	+ 3 years	+ 8 years
IRBc 4020	£105	£210
IRBc 4121	£105	£210
IRBci 4150	£125	£250
IRBci 4550/4551	£125	£250
IRBci 4571	£150	£300
IRBci 5150/5151	£125	£250
IRBci 5170/5180	£150	£300
IRBcx 5170/5171	£150	£300
IRBd 4020	£105	£210
IRBd 4050	£125	£250
IRBd 4120	£105	£210
IRBd 4150	£125	£250
IRBd 4550	£125	£250
IRBd 4570	£150	£300
IRBd 5120/5121	£105	£210
IRBd 5150/5151	£125	£250
IRBd 5170	£150	£300
IRBdi 5180/5181	£150	£300
IRBe 5120	£105	£210
IRbi 3951	£125	£250
IRBPci/cx 5170	£150	£300
IRBPdi 5170	£150	£300
IRBSd 5120/5121	£105	£210
IRBSe 5120	£105	£210
IRc(i) 3950	£125	£250
IRd 3900	£95	£190
IRd 4000	£95	£190
IRd 4021	£105	£210
IRd 4120/4121	£105	£210
IRd 4520	£105	£210
IRd 5100	£95	£190
IRDdi 5120/5121	£105	£210
IRDe 5120	£105	£210
IRe 3901	£95	£190
IRe 4100/4101	£95	£190
IRe 4520/4521	£105	£210
IRe 5100/5101	£95	£190
IRf 5101	£95	£190
IRSe 3900/3901	£95	£190
IRSe 4100/4101	£95	£190
SIbAi 3950	£125	£250

Under counter fridges

Freestanding

Model	+ 3 years	+ 8 years
Rc 1400	£60	£120
Rc 1401	£60	£120
Rci 1620	£75	£150
Rci 1621	£75	£150
Rd 1200	£60	£120
Rd 1201	£60	£120
Rd 1400	£60	£120
Rd 1400N	£60	£120
Rd 1401	£60	£120
Rdi 1620	£75	£150
Rdi 1621	£75	£150
Re 1000	£60	£120
Re 1200	£60	£120
Re 1201	£60	£120
Re 1400	£60	£120
Re 1400N	£60	£120
Re 1401	£60	£120
Rsdci 1620	£75	£150
Rsdci 1621	£75	£150
Rsve 1201	£60	£120
T 1400	£60	£120
T 1404	£60	£120
T 1410	£75	£150
T 1414	£75	£150
T 1504	£60	£120
T 1710	£75	£150
T 1714	£75	£150
T 1810	£75	£150
TP 1410	£75	£150
TP 1414	£75	£150
TP 1760	£85	£170
TP 1764	£85	£170
TPesf 1710	£85	£170
TPesf 1714	£85	£170
TX 1021	£75	£150

Built-under

Model	+ 3 years	+ 8 years
SUIB 1550	£150	£300
UIK 1510	£105	£210
UIK 1514	£105	£210
UIKo 1550	£150	£300
UIKo 1560	£150	£300
UIKP 1550	£125	£250
UIKP 1554	£125	£250
UKes 1752	£125	£250

Tall freezers

Freestanding

Model	+ 3 years	+ 8 years
FNa 6625	£105	£210
FNb 425i	£125	£250
FNb 465i	£125	£250
FNb 505i	£125	£250
FNc 467i	£150	£300
FNc 507i	£150	£300
FNc 5277	£150	£300
FNc 527i	£150	£300
FNc 6625	£105	£210
FNc 667i	£150	£300
FNc 7026	£105	£210
FNc 707i	£150	£300
FNc 7227	£105	£210
FNc 7277	£150	£300
FNc 727i	£150	£300
FNd 4254	£125	£250
FNd 4625	£105	£210
FNd 4655	£125	£250
FNd 465i	£125	£250
FNd 5026	£105	£210
FNd 5056	£125	£250
FNd 505i	£125	£250
FNd 522i	£105	£210
FNd 525i	£125	£250
FNd 7026	£105	£210
FNe 4204	£95	£190
FNe 4605	£95	£190
FNe 4625	£105	£210
FNe 5006	£95	£190
FNe 5026	£105	£210
FNe 5207	£95	£190
FNe 5227	£105	£210
FNf 4204	£95	£190
FNsdd 5297	£150	£300
FNsdd 529i	£150	£300

Built-in

Model	+ 3 years	+ 8 years
IFd 3904	£95	£190
IFe 3904	£95	£190
IFNci 3954	£125	£250
IFNd 3503	£95	£190
IFNd 3924	£105	£210
IFNd 3954	£125	£250
IFNe 3503	£95	£190
IFNe 3553	£125	£250
IFNe 3924	£105	£210
IFsd 3904	£95	£190
IFse 3904	£95	£190
SIFNAd 5188	£150	£300
SIFNAe 5188	£150	£300
SIFNci 3954	£125	£250
SIFNd 4155	£125	£250
SIFNd 4556	£125	£250
SIFNdi 4155	£125	£250
SIFNdi 4556	£125	£250
SIFNdi 5178	£150	£300
SIFNdi 5188	£150	£300
SIFNdx 5188	£150	£300
SIFNe 5108	£95	£190
SIFNe 5128	£105	£210
SIFNe 5188	£150	£300
SIFNf 5108	£95	£190
SIFNSe 5128	£105	£210
SIFNSf 5128	£105	£210

Under counter freezers

Freestanding

Model	+ 3 years	+ 8 years
Fc 1404	£60	£120
Fc 1404N	£60	£120
Fci 1624	£75	£150
Fd 1404	£60	£120
Fdi 1624	£75	£150
Fe 1003	£60	£120
Fe 1404	£60	£120
Fe 1404N	£60	£120
Fe 1414	£60	£120
FNdi 1624	£75	£150
FNe 1404	£60	£120
FNe 1404N	£60	£120
FNsddi 1624	£75	£150
Fsve 1404	£60	£120
G 1213	£60	£120
G(sl) 1223	£60	£120
GNP 1066	£85	£170
GP 1213	£60	£120
GP 1376	£75	£150
GP 1476	£75	£150
GP 1486	£85	£170
GX 823	£75	£150

Built in

Model	+ 3 years	+ 8 years
SUIG 1514	£105	£210
SUIGN 1554	£125	£250

Chest freezers

Model	+ 3 years	+ 8 years
CFe 2500	£95	£190

Wine cabinets

Freestanding

Model	+ 3 years	+ 8 years
WKb 1812	£185	£370
WKEes 553	£220	£440
WKEgb 582	£220	£440
WKes 653	£220	£440
WPbl 4201	£85	£170
WPbl 4601	£85	£170
WPbl 5001	£85	£170
WPbli 5031	£160	£320
WPbli 5231	£160	£320
WPbsi 5052	£185	£370
WPbsi 5252	£185	£370
WPgbi 5272	£220	£440
WPgbi 5273	£220	£440
WPgbi 5283	£220	£440
WPgbi 7472	£220	£440
WPgbi 7473	£220	£440
WPgbi 7483	£220	£440
WPsd 4652	£185	£370
WPsd 5252	£185	£370
WSbl 4201	£85	£170
WSbl 4601	£85	£170
WSbl 5001	£85	£170
WSbli 5031	£160	£320
WSbli 5231	£160	£320
WSbli 7731	£160	£320
WSbsi 5252	£185	£370
WTes 1672	£220	£440
WTes 5872	£220	£440
WTes 5972	£220	£440

Built in

Model	+ 3 years	+ 8 years
EWTdf 1653	£220	£440
EWTdf 2353	£220	£440
EWTdf 3553	£220	£440
EWTgb 1683	£220	£440
EWTgb 2383	£220	£440
EWTgb 3583	£220	£440
EWTgw 1683	£220	£440
EWTgw 2383	£220	£440
EWTgw 3583	£220	£440
UWKes 1752	£220	£440
UWTes 1672	£220	£440
UWTgb 1682	£220	£440
WKEes 553	£220	£440
WKEgb 582	£220	£440
WKEgw 582	£220	£440
UWKes 1752	£185	£370
UWTes 1672	£220	£440
UWTgb 1682	£220	£440

Humidor

Model	+ 3 years	+ 8 years
ZKes 453	£220	£440

Outdoor fridge

Model	+ 3 years	+ 8 years
OKes 1750	£185	£370