

IRBci 5171


Integratable refrigerator with
BioFresh Professional
△ Peak

↑
178
↓



Made
in Germany

Bio Fresh Professional	Hydro Breeze	Super Cool	Glass & Steel Interior	Light Tower	Soft System	Soft-Telescopic Rails	Super Silent	Touch & Swipe Display	Smart Device
------------------------	--------------	------------	------------------------	-------------	-------------	-----------------------	--------------	-----------------------	--------------

Type	Fixed door / Built-in appliance for integrated use
Energy efficiency class	
Energy consump. in 365 days/24 h	131 / 0.358 kWh
Total volume	276 l ¹
Volume	Fridge: 151.6 l / BioFresh compartment: 97.9 l / Freezer compartment: 27 l
Noise output / noise class	33 dB(A) / B, SuperSilent
Climate rating	SN-T (+10 °C to +43 °C)
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55
Controls	
Type of Control	2.4" TFT-colour display, touch- & swipe Digital temperature display for fridge and BioFresh compartments
Interface	SmartDeviceBox included
Door alarm	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling with Ambilight and FreshAir activated charcoal filter
Door interior	Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish
VarioSafe	yes
Bottle shelf	Integrated bottle shelf
Interior light	LED light panel, both sides LED BioFresh lighting
BioFresh-Safes	3 BioFresh safes with SoftTelescopic, 2 Fruit & Vegetable safes, of which 1 with HydroBreeze, 1 Meat & Dairy safe, 1 FlexSystem
Freezer compartment	
Star coding	****
Freezing capacity in 24 hours	2 kg
Temperature rise time	14 h
Ice cube tray	yes
Features	
Door hinges	Reversible door hinging, right
Replaceable door seal	yes
EAN-Nr.	4016803114772

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.