

# ECBN 9673

Integrated fridge-freezer with  
BioFresh and NoFrost

← 91 → 203  
↑  
↓



G

German Engineering

reddot

winner 2023

MONOLITH

Bio Fresh

Duo Cooling

No Frost

Super Cool

Super Frost

Infinity Light

Infinity Spring

Ice Maker

Soft System

Soft-Telescopic Rails

Infinity Swipe

Smart Device

|                                       |   |  |  |
|---------------------------------------|---|--|--|
| Type                                  | Fixed door / Built-in appliance for integrated use                          |  |  |
| Energy efficiency class               | <div><div>F</div><div>A</div><div>G</div></div>                             |  |  |
| Energy consump. in 365 days/24 h      | 431 / 1.180 kWh   |  |  |
| Total volume                          | 510 l <sup>1</sup>  |  |  |
| Volume                                | Fridge: 239.6 l / BioFresh compartment: 86.8 l / Freezer compartment: 184 l |  |  |
| Noise output / noise class            | 38 dB(A) / C  |  |  |
| Climate rating                        | SN-T (+10 °C to +43 °C)   |  |  |
| Aperture dimensions in cm (h / w / d) | 213.4 / 91-91.4 / min. 63.5   |  |  |
| Controls                              |   |  |  |
| Type of Control                       | 3.5" LCD colour display, behind the door, touch- & swipe                    |  |  |
|                                       | Digital temperature display   |  |  |
| Interface                             | SmartDeviceBox included   |  |  |
| Door alarm                            | acoustic  |  |  |
| Alarm in case of malfunction          | visual and acoustic   |  |  |
| Fridge compartment                    |   |  |  |
| Forced air system                     | PowerCooling with FreshAir activated charcoal filter                        |  |  |
| Door interior                         | Comfort GlassLine storage shelves with bottle holder                        |  |  |
| Interior light                        | InfinityLight   |  |  |
|                                       | LED BioFresh lighting   |  |  |
| Freezer compartment                   |   |  |  |
| Star coding                           | ****  |  |  |
| Interior light                        | LED lighting  |  |  |
| Number of drawers                     | 2 drawers on telescopic rails   |  |  |
| IceMaker                              | with fixed 3/4" water connection  |  |  |
| Ice cube capacity in 24 hours         | 1.4 kg  |  |  |
| Ice cube stock                        | 3.0 kg  |  |  |
| Freezing capacity in 24 hours         | 12 kg   |  |  |
| Temperature rise time                 | 21 h  |  |  |
| Features                              |   |  |  |
| Door hinges                           | Reversible door hinging, right  |  |  |
| EAN-Nr.                               | 9005382196052   |  |  |

<sup>1</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.