

EWT 9275

Built-in multi-temperature wine fridge



MONOLITH



213



Type	Fixed door / integrable / decoratable
Energy efficiency class	
Energy consumpt. in 365 days/24 h	180 / 0.493 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	100 ¹
Total volume	320 l ²
Noise output / noise class	38 dB(A) / C
Climate rating	SN-ST (+10 °C to +38 °C)
Frequency / voltage	50 Hz / 220-240 V~
Connection Rating	1.5 A
Door	Isolierglastür / Stahl
Aperture dimensions in cm (h / w / d)	213.4 / 61 / min. 63.5
Controls	
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe Digital temperature display for 3 wine zones
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Features	
Temperature zones	3
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	-
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	InfinityLight
Interface	SmartDeviceBox included
Shelf material	aluminium/beech
Storage shelves	9
Pull-out shelves on telescopic rails	9
Door hinges	Reversible door hinging, right
EAN-Nr.	9005382223932

¹ These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.