

# WKes 653

Wine storage fridge  
GrandCru



- 1 Zone
- Humidity Select
- Fresh Air
- UV Protect
- Vibrate Safe
- Steel Interior
- LED
- Door Lock

← 43 →



Energy efficiency class	
Energy consump. in 365 days/24 h	78 / 0.213 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	12 <sup>1</sup>
Total volume	48 l <sup>2</sup>
Noise output / noise class	37 dB(A) / C
Climate rating	SN (+10 °C to +32 °C)
Frequency / voltage	50 Hz / 220-240 V~
Connection Rating	1.0 A
Door	Glass door with stainless steel frame and SmartSteel
Side walls	stainless steel with SmartSteel
Exterior dimensions in cm (h / w / d)	61.2 / 42.5 / 47.8
<b>Controls</b>	
Type of Control	LCD display in the interior, touch electronic
	Digital temperature display
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Child proof function	yes
<b>Features</b>	
Temperature zones	1
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	1
Forced air system	Convection cooling, fresh air supply via FreshAir activated
	charcoal filter
Interior light	LED lighting permanently on, dimmable
Design	HardLine-Design
Shelf material	Beech wood fittings
Storage shelves	3
Door hinges	Fixed door hinge, right
Drawers for accessories	1
Lock	mechanisch
EAN-Nr.	9005382235492

<sup>1</sup> These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

<sup>2</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.