

CBNstc 778i

Fridge-freezer with BioFresh Professional and NoFrost

△ Peak

← 75 →




¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.



Bio Fresh Professional	Hydro Breeze	Fish & Seafood	Duo Cooling	No Frost	Super Cool	Super Frost	Vario Temp	Glass & Steel Interior	Light Tower
Soft System	Soft-Telescopic Rails	Super Silent	Touch & Swipe Display	Vario Space	Smart Device				

Energy efficiency class	
Energy consump. in 365 days/24 h	203 / 0.556 kWh ¹
Total volume	473 l ²
Volume	Fridge: 208.7 l / BioFresh compartment: 126.2 l / Freezer compartment: 138 l
Noise output / noise class	34 dB(A) / B, SuperSilent
Climate rating	SN-T (+10 °C to +43 °C)
Door / Side walls	SmartSteel / SmartSteel
Exterior dimensions in cm (h / w / d)	201.5 / 74.7 / 67.5 ¹
Controls	
Type of Control	2.4" TFT-colour display in-the-door, touch- & swipe Digital temperature display for fridge, BioFresh and freezer compartments
Interface	SmartDeviceBox included
Door alarm	visual and acoustic
Alarm in case of malfunction	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling with Ambilight and FreshAir activated charcoal filter
Door interior	Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish
VarioSafe	yes
Bottle shelf	Integrated bottle shelf
Interior light	LED light panel, both sides LED BioFresh lighting
BioFresh-Safes	2 BioFresh safes with SoftTelescopic, 1 Fruit & Vegetable safe with HydroBreeze, 1 Fish & Seafood safe, 1 FlexSystem
Freezer compartment	
Star coding	****
Interior light	LED lighting
Number of drawers	3 drawers with telescopic partial pull-out, VarioSpace
Freezing capacity in 24 hours	10 kg
Temperature rise time	20 h
Cold storage accumulators	2
Features	
Door hinges	Reversible door hinging, right
Replaceable door seal	yes
EAN-Nr.	4016803057710

LIEBHERR