

# ICBNdi 5123

Integrated fridge-freezer with  
BioFresh and NoFrost


◆ Plus

↑  
178  
↓



Made  
in Germany

Bio Fresh	Duo Cooling	No Frost	Super Cool	Super Frost	Soft System	Super Silent	Touch Display	Vario Space	Smart Device
-----------	-------------	----------	------------	-------------	-------------	--------------	---------------	-------------	--------------

Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class			
Energy consump. in 365 days/24 h	207 / 0.567 kWh		
Total volume	244 l <sup>1</sup>		
Volume	Fridge: 103.3 l / BioFresh compartment: 70.9 l / Freezer compartment: 70 l		
Noise output / noise class	34 dB(A) / B, SuperSilent		
Climate rating	SN-ST (+10 °C to +38 °C)		
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55		
Controls			
Type of Control	LC monochrome display, touch electronic Digital temperature display for fridge and freezer compartments		
Interface	SmartDeviceBox included		
Door alarm	visual and acoustic		
Alarm in case of malfunction	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling with FreshAir activated charcoal filter		
Door interior	Plastic door shelves with bottle holder, egg tray, FlexCube		
Bottle shelf	Integrated bottle shelf		
Interior light	LED ceiling lighting		
BioFresh-Safes	2 BioFresh safes with telescopic partial pull-out, 1 Fruit & Vegetable safe, 1 Meat & Dairy safe		
Freezer compartment			
Star coding	****		
Number of drawers	3 drawers, VarioSpace		
Freezing capacity in 24 hours	6 kg		
Temperature rise time	9 h		
Ice cube tray	1		
Features			
Door hinges	Reversible door hinging, right		
Replaceable door seal	yes		
EAN-Nr.	4016803115014		

<sup>1</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.