

ECBN 9671

Integrated fridge-freezer with BioFresh and NoFrost

← 91 → 203



German Engineering

reddot winner 2023

MONOLITH

Bio Fresh

No Frost

Super Cool

Super Frost

Infinity Light

Infinity Spring

Ice Maker

Soft System

Soft-Telescopic Rails

Infinity Swipe

Smart Device

Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class	<div><div></div><div>D</div><div>A</div><div>G</div></div>		
Energy consump. in 365 days/24 h	271 / 0.742 kWh		
Total volume	514 l ¹		
Volume	Fridge: 242.6 l / BioFresh compartment: 87.1 l / Freezer compartment: 184 l		
Noise output / noise class	38 dB(A) / C		
Climate rating	SN-T (+10 °C to +43 °C)		
Aperture dimensions in cm (h / w / d)	213.4 / 91-91.4 / min. 63.5		
Controls			
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe		
	Digital temperature display		
Interface	SmartDeviceBox included		
Door alarm	acoustic		
Alarm in case of malfunction	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling with FreshAir activated charcoal filter		
Door interior	Comfort GlassLine storage shelves with bottle holder		
Interior light	InfinityLight		
	LED BioFresh lighting		
Freezer compartment			
Star coding	****		
Interior light	LED lighting		
Number of drawers	2 drawers, of which 1 set on telescopic rails		
IceMaker	with fixed 3/4" water connection		
Ice cube capacity in 24 hours	1.6 kg		
Ice cube stock	3.0 kg		
Freezing capacity in 24 hours	12 kg		
Temperature rise time	21 h		
Features			
Door hinges	Reversible door hinging, right		
EAN-Nr.	9005382195994		

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.