

WTes 1672

Multi-temperature wine fridge
Vinidor



- 2 Zones
- Humidity Select
- Fresh Air
- UV Protect
- Vibrate Safe
- Steel Interior
- LED
- Telescopic Rails
- Door Lock

← 60 →



Energy efficiency class	
Energy consump. in 365 days/24 h	139 / 0.380 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	34 ¹
Total volume	95 l ²
Noise output / noise class	36 dB(A) / C
Climate rating	SN-ST (+10 °C to +38 °C)
Frequency / voltage	50 Hz / 220-240 V~
Connection Rating	1.5 A
Door	Glass door with stainless steel frame
Side walls	stainless steel
Exterior dimensions in cm (h / w / d)	82.2 / 59.8 / 57.5
Controls	
Type of Control	MagicEye control in interior, touch electronic Digital temperature display for both wine zones
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Child proof function	yes
Features	
Temperature zones	2
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	2
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	LED lighting permanently on, dimmable
Handle	Aluminium handle
Design	HardLine-Design
Shelf material	Beech wood fittings
Storage shelves	5
Pull-out shelves on telescopic rails	3
Door hinges	Reversible door hinging, right
Lock	mechanisch
EAN-Nr.	9005382235959

¹ These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.