

ICBc 5122

Integrated fridge-freezer with
BioFresh and SmartFrost


◆ Plus

↑
178
↓



Made
in Germany

Bio Fresh	Duo Cooling	Smart Frost	Super Cool	Super Frost	Soft System	Super Silent	Touch Display	Vario Space	Smart Device Ready
-----------	-------------	-------------	------------	-------------	-------------	--------------	---------------	-------------	--------------------

Type	Fixed door / Built-in appliance for integrated use
Energy efficiency class	
Energy consump. in 365 days/24 h	145 / 0.397 kWh
Total volume	255 l ¹
Volume	Fridge: 130 l / BioFresh compartment: 70.9 l / Freezer compartment: 54 l
Noise output / noise class	32 dB(A) / B, SuperSilent
Climate rating	SN-T (+10 °C to +43 °C)
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55
Controls	
Type of Control	LC monochrome display, touch electronic Digital temperature display for fridge and freezer compartments
Interface	SmartDeviceBox: optional accessory
Door alarm	visual and acoustic
Alarm in case of malfunction	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling with FreshAir activated charcoal filter
Door interior	Plastic door shelves with bottle holder, egg tray, FlexCube
Bottle shelf	Integrated bottle shelf
Interior light	LED ceiling lighting
BioFresh-Safes	2 BioFresh safes with telescopic partial pull-out, 1 Fruit & Vegetable safe, 1 Meat & Dairy safe
Freezer compartment	
Star coding	****
Number of drawers	2 drawers, VarioSpace
Freezing capacity in 24 hours	4 kg
Temperature rise time	9 h
Ice cube tray	1
Features	
Door hinges	Reversible door hinging, right
Replaceable door seal	yes
EAN-Nr.	4016803114932

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.