

# EKB 9471

Integrated refrigerator with BioFresh




MONOLITH



213



Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class			
Energy consump. in 365 days/24 h	171 / 0.468 kWh		
Total volume	440 l <sup>1</sup>		
Volume	Fridge: 250.5 l / BioFresh compartment: 190.3 l		
Noise output / noise class	37 dB(A) / C		
Climate rating	SN-T (+10 °C to +43 °C)		
Aperture dimensions in cm (h / w / d)	213.4 / 76-76.2 / min. 63.5		
Controls			
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe Digital temperature display for fridge and BioFresh compartments		
Interface	SmartDeviceBox included		
Alarm in case of malfunction	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling		
Door interior	Plastic door shelves with bottle holder		
Interior light	InfinityLight LED BioFresh lighting		
BioFresh-Safes	2 BioFresh safes with SoftTelescopic set variably to either HydroSafe or DrySafe, 1 FlexSystem		
Features			
Door hinges	Reversible door hinging, right		
EAN-Nr.	9005382233757		

<sup>1</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.