

ICBbi 5152

Integrated fridge-freezer with BioFresh and SmartFrost


Prime

↑
178
↓



Made in Germany

| | | | | | | | | | |
|-----------------------|-------------|--------------|------------|-------------|------------------------|-------------|-------------|-----------------------|--------------|
| Bio Fresh | Duo Cooling | Smart Frost | Super Cool | Super Frost | Glass & Steel Interior | Light Tower | Soft System | Soft-Telescopic Rails | Ultra Silent |
| Touch & Swipe Display | Vario Space | Smart Device | | | | | | | |

| | |
|---------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| Type | Fixed door / Built-in appliance for integrated use |
| Energy efficiency class |  |
| Energy consump. in 365 days/24 h | 116 / 0.317 kWh |
| Total volume | 256 l ¹ |
| Volume | Fridge: 131 l / BioFresh compartment: 71.2 l / Freezer compartment: 54 l |
| Noise output / noise class | 29 dB(A) / A, UltraSilent |
| Climate rating | SN-T (+10 °C to +43 °C) |
| Aperture dimensions in cm (h / w / d) | 177.2-178.8 / 56-57 / min. 55 |
| Controls | |
| Type of Control | 2.4" TFT-colour display, touch- & swipe Digital temperature display for fridge, BioFresh and freezer compartments |
| Interface | SmartDeviceBox included |
| Door alarm | visual and acoustic |
| Alarm in case of malfunction | visual and acoustic |
| Fridge compartment | |
| Forced air system | PowerCooling with FreshAir activated charcoal filter |
| Door interior | Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish |
| VarioSafe | yes |
| Bottle shelf | Integrated bottle shelf |
| Interior light | LED light panel, both sides LED BioFresh lighting |
| BioFresh-Safes | 2 BioFresh safes with SoftTelescopic, 1 Fruit & Vegetable safe, 1 Meat & Dairy safe |
| Freezer compartment | |
| Star coding | **** |
| Interior light | LED lighting |
| Number of drawers | 2 drawers, VarioSpace |
| Freezing capacity in 24 hours | 4 kg |
| Temperature rise time | 9 h |
| Cold storage accumulators | 2 |
| Ice cube tray | yes |
| Features | |
| Door hinges | Reversible door hinging, right |
| Replaceable door seal | yes |
| EAN-Nr. | 4016803115090 |

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.