

CBNbsd 578i





Fridge-freezer with BioFresh Professional and NoFrost
△ Peak

← 60 →



¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

		 reddot winner 2022								
Bio Fresh Professional	Hydro Breeze	Fish & Seafood	Duo Cooling	No Frost	Super Cool	Super Frost	Vario Temp	Glass & Steel Interior	Light Tower	
Soft System	Soft-Telescopic Rails	Super Silent	Touch & Swipe Display	Vario Space	Smart Device					
Energy efficiency class										
Energy consump. in 365 days/24 h										225 / 0.616 kWh ¹
Total volume										359 l ²
Volume										Fridge: 160.3 l / BioFresh compartment: 96.1 l / Freezer compartment: 103 l
Noise output / noise class										34 dB(A) / B, SuperSilent
Climate rating										SN-T (+10 °C to +43 °C)
Door / Side walls										BlackSteel / BlackSteel
Exterior dimensions in cm (h / w / d)										201.5 / 59.7 / 67.5 ¹
Controls										
Type of Control										2.4" TFT-colour display in-the-door, touch- & swipe Digital temperature display for fridge, BioFresh and freezer compartments
Interface										SmartDeviceBox included
Door alarm										visual and acoustic
Alarm in case of malfunction										visual and acoustic
Fridge compartment										
Forced air system										PowerCooling with Ambilight and FreshAir activated charcoal filter
Door interior										Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish
VarioSafe										yes
Bottle shelf										Integrated bottle shelf
Interior light										LED light panel, both sides LED BioFresh lighting
BioFresh-Safes										2 BioFresh safes with SoftTelescopic, 1 Fruit & Vegetable safe with HydroBreeze, 1 Fish & Seafood safe, 1 FlexSystem
Freezer compartment										
Star coding										****
Interior light										LED lighting
Number of drawers										3 drawers, of which 2 with telescopic partial pull-out, VarioSpace
Freezing capacity in 24 hours										10 kg
Temperature rise time										20 h
Cold storage accumulators										2
Features										
Door hinges										Reversible door hinging, right
Replaceable door seal										yes
EAN-Nr.										4016803057475

LIEBHERR