

IRBd 5121

Integrated refrigerator with BioFresh


Plus

178



Made in Germany

- Bio Fresh
- Super Cool
- Soft System
- Super Silent
- Touch Display
- Smart Device Ready

Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class			
Energy consump. in 365 days/24 h	164 / 0.449 kWh		
Total volume	275 l ¹		
Volume	Fridge: 150.2 l / BioFresh compartment: 98.2 l / Freezer compartment: 27 l		
Noise output / noise class	34 dB(A) / B, SuperSilent		
Climate rating	SN-T (+10 °C to +43 °C)		
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55		
Controls			
Type of Control	LC monochrome display, touch electronic Digital temperature display for fridge compartment		
Interface	SmartDeviceBox: optional accessory		
Door alarm	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling with FreshAir activated charcoal filter		
Door interior	Plastic door shelves with bottle holder, egg tray, FlexCube		
Bottle shelf	Integrated bottle shelf		
Interior light	LED light stripe, both sides		
BioFresh-Safes	3 BioFresh safes with telescopic partial pull-out, 2 Fruit & Vegetable safes, 1 Meat & Dairy safe		
Freezer compartment			
Star coding	****		
Freezing capacity in 24 hours	2 kg		
Temperature rise time	14 h		
Ice cube tray	1		
Features			
Door hinges	Reversible door hinging, right		
Replaceable door seal	yes		
EAN-Nr.	4016803113638		

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.