

EKB 9471

Integrated refrigerator with BioFresh




MONOLITH



213



Type	Built-in appliance for integrated use	
Energy efficiency class		
Energy consump. in 365 days/24 h	171 / 0.468 kWh	
Total volume	440 l ¹	
Volume	Fridge: 250.5 l / BioFresh compartment: 190.3 l	
Noise output / noise class	37 dB(A) / C	
Climate rating	SN-T (+10 °C to +43 °C)	
Aperture dimensions in cm (h / w / d)	213.4 / 76.2 / min. 63.5	
Controls		
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe Digital temperature display for fridge and BioFresh compartments	
Interface	SmartDeviceBox included	
Fridge compartment		
Door interior	Plastic door shelves with bottle holder	
Interior light	InfinityLight LED BioFresh lighting	
BioFresh-Safes	2 BioFresh safes with SoftTelescopic set variably to either HydroSafe or DrySafe, 1 FlexSystem	
Features		
EAN-Nr.	9005382274651	

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.