

CBNsdA 575i

Fridge-freezer with BioFresh and NoFrost

Prime

← 60 →



¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

<div><div><div>G</div><div>German Engineering</div></div></div>		<div><div>BLU</div><div>Performance</div></div>	
<div><div>Bio Fresh</div><div>Duo Cooling</div><div>No Frost</div><div>Super Cool</div><div>Super Frost</div><div>Vario Temp</div><div>Glass & Steel Interior</div><div>Light Tower</div><div>Easy Twist Ice</div><div>Soft System</div><div>Soft-Telescopic Rails</div><div>Super Silent</div><div>Touch & Swipe Display</div><div>Vario Space</div><div>Smart Device</div></div>			
Energy efficiency class		<div><div>A</div><div>A</div><div>G</div></div>	
Energy consump. in 365 days/24 h		116 / 0.317 kWh ¹	
Total volume		362 l ²	
Volume		Fridge: 159.9 l / BioFresh compartment: 99.2 l / Freezer compartment: 103 l	
Noise output / noise class		33 dB(A) / B, SuperSilent	
Climate rating		SN-T (+10 °C to +43 °C)	
Door / Side walls		SmartSteel / silver	
Exterior dimensions in cm (h / w / d)		201.5 / 59.7 / 67.5 ¹	
Controls			
Type of Control		2.4" TFT-colour display behind-the-door, touch- & swipe Digital temperature display for fridge, BioFresh and freezer compartments	
Interface		SmartDeviceBox included	
Door alarm		visual and acoustic	
Alarm in case of malfunction		visual and acoustic	
Fridge compartment			
Forced air system		PowerCooling with FreshAir activated charcoal filter	
Door interior		Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish	
VarioSafe		yes	
Bottle shelf		Integrated bottle shelf	
Interior light		LED light panel, both sides LED BioFresh lighting	
BioFresh-Safes		2 BioFresh safes with SoftTelescopic, 1 Fruit & Vegetable safe, 1 Meat & Dairy safe	
Freezer compartment			
Star coding		****	
Interior light		LED lighting	
Number of drawers		3 drawers, VarioSpace	
Freezing capacity in 24 hours		10 kg	
Temperature rise time		20 h	
Cold storage accumulators		2	
Features			
Door hinges		Reversible door hinging, right	
Replaceable door seal		yes	
EAN-Nr.		4016803129370	

LIEBHERR