

# EWTdf 2353

Built-in multi-temperature wine fridge  
Vinidor

↑  
122  
↓



2 Zones	Humidity Select	Fresh Air	UV Protect	Vibrate Safe	Steel Interior	LED	Soft System	Telescopic Rails	Super Silent
Presentation Board									

Type	Fixed door / Built-in appliance with decor panel
Energy efficiency class	
Energy consump. in 365 days/24 h	151 / 0.413 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	48 <sup>1</sup>
Total volume	159 l <sup>2</sup>
Noise output / noise class	32 dB(A) / B, SuperSilent
Climate rating	SN-ST (+10 °C to +38 °C)
Frequency / voltage	50/60 Hz / 220-240 V~
Connection Rating	1.5 A
Door	Insulated glass door prepared for cabinet door or stainless steel front
Aperture dimensions in cm (h / w / d)	122-123.6 / 56-57 / min. 55
Controls	
Type of Control	LCD display, touch electronic Digital temperature display for both wine zones
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Child proof function	yes
Features	
Temperature zones	2
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	2
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	LED lighting permanently on, dimmable
Shelf material	Beech wood fittings
Storage shelves	6
Pull-out shelves on telescopic rails	4
Door hinges	Reversible door hinging, right
EAN-Nr.	9005382278932

<sup>1</sup> These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

<sup>2</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.