

RBsfd 5221

Refrigerator with BioFresh

◆ Plus



← 60 →



Energy efficiency class	
Energy consump. in 365 days/24 h	164 / 0.449 kWh ¹
Total volume	351 l ²
Volume	Fridge: 223 l / BioFresh compartment: 93.9 l / Freezer compartment: 34 l
Noise output / noise class	35 dB(A) / B, SuperSilent
Climate rating	SN-T (+10 °C to +43 °C)
Door / Side walls	SteelFinish / silver
Exterior dimensions in cm (h / w / d)	185.5 / 59.7 / 67.5 ¹
Controls	
Type of Control	LC monochrome display behind-the-door, touch electronic Digital temperature display for fridge compartment
Interface	SmartDeviceBox: optional accessory
Door alarm	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling with FreshAir activated charcoal filter
Door interior	Plastic door shelves with bottle holder, egg tray, FlexCube
Bottle shelf	Integrated bottle shelf
Interior light	LED light one-sided
BioFresh-Safes	2 BioFresh safes with telescopic partial pull-out, 1 Fruit & Vegetable safe, 1 Meat & Dairy safe
Freezer compartment	
Star coding	****
Freezing capacity in 24 hours	5 kg
Temperature rise time	16 h
Features	
Door hinges	Reversible door hinging, right
Replaceable door seal	yes
EAN-Nr.	4016803130611

¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.