

ECBN 9673

Integrated fridge-freezer with BioFresh and NoFrost

← 91 → 203



German Engineering

reddot winner 2023

MONOLITH

Bio Fresh

Duo Cooling

No Frost

Super Cool

Super Frost

Infinity Light

Infinity Spring

Ice Maker

Soft System

Infinity Swipe

Smart Device

Type	Fixed door / Built-in appliance for integrated use	
Energy efficiency class	<div><div>F</div><div>A</div><div>G</div></div>	
Energy consump. in 365 days/24 h	431 / 1.180 kWh	
Total volume	510 l ¹	
Volume	Fridge: 239.6 l / BioFresh compartment: 86.8 l / Freezer compartment: 184 l	
Noise output / noise class	38 dB(A) / C	
Climate rating	SN-T (+10 °C to +43 °C)	
Aperture dimensions in cm (h / w / d)	213.4 / 91.4 / min. 63.5	
Controls		
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe	
	Digital temperature display	
Interface	SmartDeviceBox included	
Door alarm	acoustic	
Alarm in case of malfunction	visual and acoustic	
Fridge compartment		
Forced air system	PowerCooling with FreshAir activated charcoal filter	
Door interior	Comfort GlassLine storage shelves with bottle holder	
Interior light	InfinityLight	
	LED BioFresh lighting	
Freezer compartment		
Star coding	****	
Interior light	LED lighting	
Number of drawers	2 drawers on telescopic rails	
IceMaker	with fixed 3/4" water connection	
Ice cube capacity in 24 hours	1.4 kg	
Ice cube stock	3.0 kg	
Freezing capacity in 24 hours	12 kg	
Temperature rise time	21 h	
Features		
Door hinges	Reversible door hinging, right	
EAN-Nr.	9005382274576	

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

Subject to modification. Timestamp: 14.05.2024

LIEBHERR