

ICBNbsci 5173

Integrated fridge-freezer with
BioFresh Professional and NoFrost
△ Peak

↑
178
↓



¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.



Bio Fresh Professional	Hydro Breeze	Duo Cooling	No Frost	Super Cool	Super Frost	Vario Temp	Glass & Steel Interior	Light Tower	Ice Maker
Soft System	Soft-Telescopic Rails	Super Silent	Touch & Swipe Display	Vario Space	Smart Device				

Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class	<div><div>C</div><div>A + G</div></div>		
Energy consump. in 365 days/24 h	166 / 0.454 kWh		
Total volume	246 l ¹		
Volume	Fridge: 104.5 l / BioFresh compartment: 70.7 l / Freezer compartment: 71 l		
Noise output / noise class	33 dB(A) / B, SuperSilent		
Climate rating	SN-ST (+10 °C to +38 °C)		
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55		
Controls			
Type of Control	2.4" TFT-colour display, touch- & swipe Digital temperature display for fridge, BioFresh and freezer compartments		
Interface	SmartDeviceBox included		
Door alarm	visual and acoustic		
Alarm in case of malfunction	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling with Ambilight and FreshAir activated charcoal filter		
Door interior	Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish		
Bottle shelf	Integrated bottle shelf		
Interior light	LED light panel, both sides LED BioFresh lighting		
BioFresh-Safes	2 BioFresh safes with SoftTelescopic, 1 Fruit & Vegetable safe with HydroBreeze, 1 Meat & Dairy safe, 1 FlexSystem		
Freezer compartment			
Star coding	****		
Interior light	LED lighting		
Number of drawers	3 drawers, of which 2 with telescopic partial pull-out, VarioSpace		
IceMaker	with fixed 3/4" water connection, LED lighting		
Ice cube capacity in 24 hours	1.2 kg		
Ice cube stock	5.0 kg		
Freezing capacity in 24 hours	6 kg		
Temperature rise time	9 h		
Cold storage accumulators	2		
Features			
Door hinges	Reversible door hinging, right		
Replaceable door seal	yes		
EAN-Nr.	4016803120971		