

EWT 9175

Built-in multi-temperature wine fridge



MONOLITH



213



¹ These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

Type	Fixed door / integrable / decoratable
Energy efficiency class	
Energy consumpt. in 365 days/24 h	154 / 0.421 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	75 ¹
Total volume	232 l ²
Noise output / noise class	38 dB(A) / C
Climate rating	SN-ST (+10 °C to +38 °C)
Frequency / voltage	50 Hz / 220-240 V~
Connection Rating	1.5 A
Door	Isolierglastür / Stahl
Aperture dimensions in cm (h / w / d)	213.4 / 45.7 / min. 63.5
Controls	
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe Digital temperature display for both wine zones
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Features	
Temperature zones	2
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	-
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	InfinityLight
Interface	SmartDeviceBox included
Shelf material	aluminium/beech
Storage shelves	9
Pull-out shelves on telescopic rails	9
Door hinges	Reversible door hinging, right
EAN-Nr.	9005382196090

LIEBHERR