

IRBd 5120

Integrated refrigerator with BioFresh

◆ Plus



↑
178
↓



Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class			
Energy consump. in 365 days/24 h	125 / 0.342 kWh		
Total volume	294 l ¹		
Volume	Fridge: 195.8 l / BioFresh compartment: 98.2 l		
Noise output / noise class	33 dB(A) / B, SuperSilent		
Climate rating	SN-T (+10 °C to +43 °C)		
Aperture dimensions in cm (h / w / d)	177.2-178.8 / 56-57 / min. 55		
Controls			
Type of Control	LC monochrome display, touch electronic Digital temperature display for fridge compartment		
Interface	SmartDeviceBox: optional accessory		
Door alarm	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling with FreshAir activated charcoal filter		
Door interior	Plastic door shelves with bottle holder, egg tray, FlexCube		
Bottle shelf	Integrated bottle shelf		
Interior light	LED ceiling lighting		
BioFresh-Safes	3 BioFresh safes with telescopic partial pull-out, 2 Fruit & Vegetable safes, 1 Meat & Dairy safe		
Features			
Door hinges	Reversible door hinging, right		
Replaceable door seal	yes		
EAN-Nr.	4016803113591		

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.