

# CBNste 8872

Fridge-freezer with BioFresh and NoFrost



← 91 →



Energy efficiency class	
Energy consump. in 365 days/24 h	320 / 0.876 kWh
Total volume	522 l <sup>1</sup>
Volume	Fridge: 292.2 l / BioFresh compartment: 70.6 l / Freezer compartment: 160 l
Noise output / noise class	42 dB(A) / D
Climate rating	SN-T (+10 °C to +43 °C)
Door / Side walls	Stainless steel / stainless steel
Exterior dimensions in cm (h / w / d)	202.7 / 91 / 61.5
Controls	
Type of Control	MagicEye behind-the-door, touch electronic Digital temperature display for fridge and freezer compartments
Door alarm	acoustic
Alarm in case of malfunction	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling
Door interior	Premium GlassLine storage shelves with bottle holder
Bottle shelf	Integrated bottle shelf
Interior light	LED light column on both sides and LED ceiling lighting LED BioFresh lighting
BioFresh-Safes	2 BioFresh safes on telescopic rails set variably to either HydroSafe or DrySafe
Freezer compartment	
Star coding	****
Interior light	LED lighting
Number of drawers	2 drawers on telescopic rails
IceMaker	with fixed 3/4" water connection
Ice cube capacity in 24 hours	1.4 kg
Ice cube stock	2.7 kg
Freezing capacity in 24 hours	9.5 kg
Temperature rise time	16.9 h
Features	
Door hinges	Door hinge left fixed/right fixed
Replaceable door seal	yes
EAN-Nr.	9005382270875

<sup>1</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.