

CBNbsa 575i

Fridge-freezer with BioFresh and NoFrost

Prime

← 60 →



¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

<div><div><div>G</div><div>German Engineering</div></div></div>														<div><div>BLU</div><div>Performance</div></div>																													
Bio Fresh				Duo Cooling				No Frost				Super Cool				Super Frost				Vario Temp				Glass & Steel Interior				Light Tower				Easy Twist Ice				Soft System							
Soft-Telescopic Rails				Super Silent				Touch & Swipe Display				Vario Space				Smart Device																											
Energy efficiency class																																<div><div>A</div><div>A</div><div>G</div></div>											
Energy consump. in 365 days/24 h																																116 / 0.317 kWh ¹											
Total volume																																362 l ²											
Volume																																Fridge: 159.9 l / BioFresh compartment: 99.2 l / Freezer compartment: 103 l											
Noise output / noise class																																33 dB(A) / B, SuperSilent											
Climate rating																																SN-T (+10 °C to +43 °C)											
Door / Side walls																																BlackSteel / BlackSteel											
Exterior dimensions in cm (h / w / d)																																201.5 / 59.7 / 67.5 ¹											
Controls																																											
Type of Control																																2.4" TFT-colour display behind-the-door, touch- & swipe				Digital temperature display for fridge, BioFresh and freezer compartments							
Interface																																SmartDeviceBox included											
Door alarm																																visual and acoustic											
Alarm in case of malfunction																																visual and acoustic											
Fridge compartment																																											
Forced air system																																PowerCooling with FreshAir activated charcoal filter											
Door interior																																Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish											
VarioSafe																																yes											
Bottle shelf																																Integrated bottle shelf											
Interior light																																LED light panel, both sides				LED BioFresh lighting							
BioFresh-Safes																																2 BioFresh safes with SoftTelescopic,				1 Fruit & Vegetable safe,				1 Meat & Dairy safe			
Freezer compartment																																											
Star coding																																****											
Interior light																																LED lighting											
Number of drawers																																3 drawers, VarioSpace											
Freezing capacity in 24 hours																																10 kg											
Temperature rise time																																20 h											
Cold storage accumulators																																2											
Features																																											
Door hinges																																Reversible door hinging, right											
Replaceable door seal																																yes											
EAN-Nr.																																4016803129356											

LIEBHERR