

CBNbsd 578i

Fridge-freezer with BioFresh Professional and NoFrost
△ Peak

← 60 →



¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

German Engineering

reddot winner 2022

BLU

Performance

Bio Fresh Professional

Hydro Breeze

Fish & Seafood

Duo Cooling

No Frost

Super Cool

Super Frost

Vario Temp

Glass & Steel Interior

Light Tower

Soft System

Soft-Telescopic Rails

Super Silent

Touch & Swipe Display

Vario Space

Smart Device

Energy efficiency class

D

A

G

Energy consump. in 365 days/24 h

225 / 0.616 kWh¹

Total volume

359 l²

Volume

Fridge: 160.3 l / BioFresh compartment: 96.1 l / Freezer compartment: 103 l

Noise output / noise class

34 dB(A) / B, SuperSilent

Climate rating

SN-T (+10 °C to +43 °C)

Door / Side walls

BlackSteel / BlackSteel

Exterior dimensions in cm (h / w / d)

201.5 / 59.7 / 67.5¹

Controls

Type of Control

2.4" TFT-colour display in-the-door, touch- & swipe
Digital temperature display for fridge, BioFresh and freezer compartments

Interface

SmartDeviceBox included

Door alarm

visual and acoustic

Alarm in case of malfunction

visual and acoustic

Fridge compartment

Forced air system

PowerCooling with Ambilight and FreshAir activated charcoal filter

Door interior

Premium GlassLine storage shelves with bottle holder, removable VarioBoxes, adjustable egg tray, butter dish

VarioSafe

yes

Bottle shelf

Integrated bottle shelf

Interior light

LED light panel, both sides
LED BioFresh lighting

BioFresh-Safes

2 BioFresh safes with SoftTelescopic,
1 Fruit & Vegetable safe with HydroBreeze,
1 Fish & Seafood safe, 1 FlexSystem

Freezer compartment

Star coding

Interior light

LED lighting

Number of drawers

3 drawers, of which 2 with telescopic partial pull-out, VarioSpace

Freezing capacity in 24 hours

10 kg

Temperature rise time

20 h

Cold storage accumulators

2

Features

Door hinges

Reversible door hinging, right

Replaceable door seal

yes

EAN-Nr.

4016803129394

LIEBHERR