

# EWTdf 2353

Built-in multi-temperature wine fridge  
Vinidor

↑  
122  
↓



2 Zones	Humidity Select	Fresh Air	UV Protect	Vibrate Safe	Steel Interior	LED	Soft System	Telescopic Rails	Super Silent
Presentation Board									

Type	Fixed door / Built-in appliance with decor panel
Energy efficiency class	
Energy consump. in 365 days/24 h	158 / 0.432 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	48 <sup>1</sup>
Total volume	158 l <sup>2</sup>
Noise output / noise class	32 dB(A) / B, SuperSilent
Climate rating	SN-ST (+10 °C to +38 °C)
Frequency / voltage	50 Hz / 220-240 V~
Connection Rating	1.0 A
Door	Insulated glass door prepared for cabinet door or stainless steel front
Aperture dimensions in cm (h / w / d)	122-123.6 / 56-57 / min. 55
<b>Controls</b>	
Type of Control	LCD display, touch electronic Digital temperature display for both wine zones
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Child proof function	yes
<b>Features</b>	
Temperature zones	2
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	2
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	LED lighting permanently on, dimmable
Shelf material	Beech wood fittings
Storage shelves	6
Pull-out shelves on telescopic rails	4
Door hinges	Reversible door hinging, right
EAN-Nr.	9005382234013

<sup>1</sup> These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

<sup>2</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.