

# ECBNe 8870

Integrated fridge-freezer with  
BioFresh and NoFrost

← 91 → 203  
↑  
↓



Type	Fixed door / Built-in appliance for integrated use		
Energy efficiency class			
Energy consump. in 365 days/24 h	320 / 0.876 kWh		
Total volume	522 l <sup>1</sup>		
Volume	Fridge: 292.2 l / BioFresh compartment: 70.6 l / Freezer compartment: 160 l		
Noise output / noise class	39 dB(A) / C		
Climate rating	SN-T (+10 °C to +43 °C)		
Aperture dimensions in cm (h / w / d)	203.2-207.2 / 91.5 / min. 63.5		
Controls			
Type of Control	MagicEye, touch electronic Digital temperature display for fridge and freezer compartments		
Door alarm	acoustic		
Alarm in case of malfunction	visual and acoustic		
Fridge compartment			
Forced air system	PowerCooling		
Door interior	Comfort GlassLine storage shelves with bottle holder		
Bottle shelf	Integrated bottle shelf		
Interior light	LED light column on both sides and LED ceiling lighting LED BioFresh lighting		
BioFresh-Safes	2 BioFresh safes on telescopic rails set variably to either HydroSafe or DrySafe		
Freezer compartment			
Star coding	****		
Interior light	LED lighting		
Number of drawers	2 drawers on telescopic rails		
IceMaker	with fixed 3/4" water connection		
Ice cube capacity in 24 hours	1.4 kg		
Ice cube stock	2.7 kg		
Freezing capacity in 24 hours	9.5 kg		
Temperature rise time	16.9 h		
Features			
Door hinges	Reversible door hinging, right		
Replaceable door seal	yes		
EAN-Nr.	9005382270912		

<sup>1</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.