

# EWTgb 3583

Built-in multi-temperature wine fridge  
Vinidor

↑  
178  
↓



2 Zones	Humidity Select	Fresh Air	UV Protect	Vibrate Safe	Steel Interior	LED	Soft System	Telescopic Rails	Tip Open
Super Silent	Presentation Board								

Type	flush / Built-in appliance for integrated use
Energy efficiency class	
Energy consump. in 365 days/24 h	165 / 0.452 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	83 <sup>1</sup>
Total volume	276 l <sup>2</sup>
Noise output / noise class	32 dB(A) / B, SuperSilent
Climate rating	SN (+10 °C to +32 °C)
Frequency / voltage	50/60 Hz / 220-240 V~
Connection Rating	1.5 A
Door	black insulated glass door, handle and frame free
Aperture dimensions in cm (h / w / d)	177.2-178 / 56-57 / min. 55
<b>Controls</b>	
Type of Control	LCD display, touch electronic Digital temperature display for both wine zones
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Child proof function	yes
<b>Features</b>	
Temperature zones	2
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	2
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	LED lighting permanently on, dimmable
Shelf material	Beech wood fittings
Storage shelves	10
Pull-out shelves on telescopic rails	8
Door hinges	Reversible door hinging, right
EAN-Nr.	9005382279038

<sup>1</sup> These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

<sup>2</sup> In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.