

EWT 9175

Built-in multi-temperature wine fridge



MONOLITH



213



Type	Fixed door / integrable / decoratable	
Energy efficiency class		
Energy consump. in 365 days/24 h	154 / 0.421 kWh	
Max. No. of Bordeaux bottles 0.75 Ltr	75 ¹	
Total volume	232 l ²	
Noise output / noise class	38 dB(A) / C	
Climate rating	SN-ST (+10 °C to +38 °C)	
Frequency / voltage	50 Hz / 220-240 V~	
Connection Rating	1.5 A	
Door	Isolierglastür / Stahl	
Aperture dimensions in cm (h / w / d)	213.4 / 45.7 / min. 63.5	
Controls		
Type of Control	3.5" LCD colour display, behind the door, touch- & swipe Digital temperature display for both wine zones	
Alarm in case of malfunction	visual and acoustic	
Door alarm	acoustic	
Features		
Temperature zones	2	
Temperature range	+5 °C to +20 °C	
Adjustable cooling circuits	-	
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter	
Interior light	InfinityLight	
Interface	SmartDeviceBox included	
Shelf material	aluminium/beech	
Storage shelves	9	
Door hinges	Reversible door hinge, left	
EAN-Nr.	9005382274118	

¹ These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

² In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.