
Claim £100 CASHBACK

Buy selected fridge-freezers
before 31st December 2024
and enjoy £100 CASHBACK.

LIEBHERR

Refrigeration and Freezing



Claim **£100 CASHBACK**

Models included



£100

Cash Back Value



CBNa 572i Plus

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 59.7 / 67.5

Total volume ¹: 360 litre



£100

Cash Back Value



CBNbda 572i Plus

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 59.7 / 67.5

Total volume ¹: 360 litre



£100

Cash Back Value



CBNbdc 778i Peak

Fridge-freezer with BioFresh Professional and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 74.7 / 67.5

Total volume ¹: 473 litre



£100

Cash Back Value



CBNc 5723 Plus

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 59.7 / 67.5

Total volume ¹: 360 litre



£100

Cash Back Value



CBNsdb 575i Prime

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 59.7 / 67.5

Total volume ¹: 362 litre



£100

Cash Back Value



CBNsdc 573i Plus

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 59.7 / 67.5

Total volume ¹: 360 litre



£100

Cash Back Value



CBNsdc 765i Prime

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 185.5 / 74.7 / 67.5

Total volume ¹: 423 litre



£100

Cash Back Value



CBNsfc 5223 Plus

Fridge-freezer with BioFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 185.5 / 59.7 / 67.5

Total volume ¹: 320 litre



£100

Cash Back Value



CNbdb 5223 Plus

Fridge-freezer with EasyFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 185.5 / 59.7 / 67.5

Total volume ¹: 330 litre



£100

Cash Back Value



CNsdb 5723 Plus

Fridge-freezer with EasyFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 201.5 / 59.7 / 67.5

Total volume ¹: 371 litre



£100

Cash Back Value



CNsdc 5203 Pure

Fridge-freezer with EasyFresh and NoFrost

Appliance Dimensions:
HxWxD (cm) 185.5 / 74.7 / 67.5

Total volume ¹: 330 litre

¹ In accordance with Regulation EU 2019/2016 we indicate the total volume to the nearest integer (rounded off) and the volume of the frozen compartments and refrigeration compartments to one place after the decimal point.

* at participating retailers only. From 01/10/2024 until 31/12/2024.

How to claim

- Claims must be made at least 30 days after purchase, but no later than 28/02/2025
- Please visit www.liebherrpromotions.co.uk to create a profile and start your claim, or check/update an existing claim. You will need: An electronic copy of your original sales receipt showing the qualifying model purchased and the date of sale, details of the retailer you purchased from, the appliance Serial Number+, your full address, and an email address which will be used to communicate claim updates to you.
- Your reward of up to £100 will be paid via a prepaid Liebherr branded MasterCard and will be approved within 30 days of receiving your completed online claim form. There is no cash alternative for this promotion.
- Promotion valid for all purchases made between 01/10/2024 and 31/12/2024, at participating retailers only, on the below products.
- All claims are subject to the promotion's Terms and Conditions.



+ Our 9-digit Serial Number can be found on the Rating Plate of your appliance. This is usually located at the bottom left, inside the appliance. To view it, you might have to remove a drawer or shelf first.